

Commentary

Open Access



Combining Cre-LoxP and single-cell sequencing technologies: insights into the extracellular vesicle cargo transfer

Michael W. Pfaffl 

Animal Physiology and Immunology, School of Life Sciences, Technical University of Munich, Freising 85354, Germany.

Correspondence to: Prof. Michael W. Pfaffl, Animal Physiology and Immunology, School of Life Sciences, Technical University of Munich, Liesel-Beckmann-Straße 1, Freising 85354, Germany. E-mail: michael.pfaffl@tum.de

How to cite this article: Pfaffl MW. Combining Cre-LoxP and single-cell sequencing technologies: insights into the extracellular vesicle cargo transfer. *Extracell Vesicles Circ Nucleic Acids* 2024;5:614-7. <https://dx.doi.org/10.20517/evcna.2024.58>

Received: 2 Aug 2024 **First Decision:** 27 Sep 2024 **Revised:** 22 Oct 2024 **Accepted:** 28 Oct 2024 **Published:** 9 Nov 2024

Academic Editors: Shenglin Huang, Yoke Peng Loh **Copy Editor:** Ting-Ting Hu **Production Editor:** Ting-Ting Hu

Abstract

The recent study from the Pogge von Strandmann group published in *Cellular and Molecular Immunology*, by Alashkar Alhamwe *et al.*, combined for the first time the Cre-LoxP recombination system with single-cell sequencing. The group monitored the tumor-derived extracellular vesicle (EV) uptake and the EV functions in the recipient non-malignant cells in a pancreatic ductal adenocarcinoma mouse model. Recombination events and EV uptake, together with resulting gene expression changes in macrophages, neutrophils, and mast cells, were detected by single-cell sequencing technology of the tumor tissue. This new approach is highly specific, as it can identify single EV recipient cells without interfering with the EV biogenesis or the phenotype.

Keywords: Extracellular vesicles, intracellular cargo transport, pancreatic ductal adenocarcinoma, tumor microenvironments, Cre-LoxP, single-cell sequencing, Bag6

TEXT

Extracellular vesicles (EVs) are diverse, nanosized, double membrane-bound structures secreted by nearly all cell types under healthy and pathological conditions. They play a crucial role in mediating intercellular communication by transferring and delivering biologically active molecules, such as various types of nucleic acids, proteins, metabolites, and lipids, to adjacent tissues or distant recipient cells, thereby eliciting specific



© The Author(s) 2024. **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License (<https://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, sharing, adaptation, distribution and reproduction in any medium or format, for any purpose, even commercially, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made.



biological responses^[1].

Despite significant advancements in EV research during the last two decades, a major challenge impeding their diagnostic and therapeutic application is the limited understanding of their *in vivo* biological functions. Determining the EV biodistribution in tissues and organs and the underlying EV kinetics has been the primary focus of most *in vivo* studies, with a possible emphasis on their role in immune recipient cells. Based on the minimal information for studies of extracellular vesicles (MISEV) 2023 guideline, fluorescence and bioluminescent tags have been used to monitor tumors and EVs *in vivo* by labeling and tracing cells^[1-4]. However, several disadvantages arise, including issues with EV labeling efficiency, deep tissue penetration, and signal-to-noise ratio. Furthermore, the type of labeling tag itself, e.g., using green fluorescent protein (GFP), can alter the EV biogenesis, the phenotype, and finally, the function. In addition, the use of additional genetic manipulation approaches, such as incorporating a tag fusion construct (e.g., CD63-GFP), may lead to increased off-target effects in non-malignant immune cells.

Today, the Cre-LoxP recombination system is widely established as a highly effective method to study EV uptake and function *in vivo*^[5-7]. To visualize the uptake of EVs by recipient cells, a reporter mouse strain such as B6.129(Cg)-Gt(ROSA)26Sortm4(ACTB-tdTomato,-eGFP)Luo/J is used. If vesicles equipped with the Cre recombinase are taken up by target cells, the excision of a LoxP-flanked stop cassette upstream of eGFP occurs and enables eGFP expression. Green fluorescent target cells can then be identified, e.g., via immunohistochemistry of tissues or upon isolation using flow cytometry^[5-7]. One recent tool to study the *in vivo* transfer of EVs is the Exomap1 transgenic mouse, which, in response to Cre recombinase, expresses an exosomal marker protein^[8]. Using *in vivo* tumor models offers a deeper understanding of the direct transfer of EV cargo between malignant cells and, furthermore, from malignant to non-malignant cells within the tumor microenvironment (TME). However, so far, the identification and characterization of the fluorescent target cells remain technically challenging.

A recent study from the Pogge von Strandmann lab combined, for the first time, the Cre-LoxP system with a state-of-the-art single-cell sequencing approach to tackle this limitation. The aim of the study was to investigate the role of the immune cell regulator and chaperone BAG6 (Bcl2-associated-athanogene 6) in EV cargo loading and function within the pancreatic cancer TME^[9]. To this end, the authors used a preclinical mouse model for pancreatic ductal adenocarcinoma (PDAC) in a BAG6 pro- or deficient background employing Cre-LoxP reporter mice. The uptake of Cre recombinase mRNA, delivered through EVs from transplanted tumor cells, was tracked at the molecular level by single-cell RNA sequencing (scRNAseq) of the tumor tissue. The beauty of the system is that the functional impact of EV uptake on the target cells can be directly analyzed at the single-cell level. This is possible by comparing the gene expression pattern of cells with and without vesicle uptake, which is indicated by recombination and, therefore, eGFP expression^[9]. Combining the Cre-LoxP system with scRNA sequencing thus allows for determining the gene profile of target cells upon EV uptake in a direct comparison with non-recombined cells in the tissue. Thus, specific EV target cell populations can be identified, as the proteins they present offer a reliable tool for investigating the functional role of EVs and their impact on cellular signaling and downstream effects on various biological processes. The validity of the method was verified through immunofluorescence staining of GFP protein in the tumor tissues of animals transplanted with Cre+ or Cre- tumor cells.

Using this approach, the authors elegantly showed that the *in vivo* EV uptake induced changes in the cellular composition of the TME in a Bag6-dependent manner. Recombination and vesicle uptake were detectable, e.g., in macrophages and neutrophils, but of note, predominantly and specifically in the mast cell compartment of Bag6-deficient tumors. Further analysis revealed that this activation was triggered by the

interaction between the Interleukin-33 receptor and Interleukin-33 presenting tumor-derived EVs. This EV-initiated cascade induced changes within the TME, generating a tumor-promoting milieu and accelerated tumor growth. *In vivo* tracking of EVs involves labeling EVs or tumor cells with fluorescent dyes or incorporating tag fusion constructs, which can affect the biogenesis, integrity, surface properties, and biological activity of EVs. To address these limitations, the Cre-LoxP-based technique offers a precise method for monitoring EV transfer in both *in vitro* and *in vivo* models via fluorescent labeling to specifically determine EV target cells^[4]. By using this approach together with single-cell sequencing, it is possible to track EV uptake *in vivo* and, at the same time, to investigate the impact of vesicle uptake on gene expression in the recipient cell. The authors applied the Pan02 model for their studies, which does not harbor the pancreatic cancer driver mutations Kirsten rat sarcoma viral oncogene homolog (KRAS), tumor protein 53 (TP53), and cyclin-dependent kinase inhibitor 2A (CDKN2A), but a SMAD family member 4 (SMAD4) mutation^[10,11]. Pan02 cells thus express wild-type p53, which is known to be involved in Bag6-dependent EV biogenesis^[12]. Therefore, it would be interesting to compare the EV phenotype of Pan02 tumors with other models with different genetic backgrounds with a view on PDAC heterogeneity. This method is not limited to PDAC mouse models but can be applied generally in research addressing the function of EVs in physiological conditions and diseases. The approach is highly specific, may identify single EV recipient cells, and unravels how tumor-derived EVs alter the phenotype of non-malignant immune cells in the TME. In the future, further developments are needed to elucidate the influence of proteins presented on the corona of the EVs, such as Interleukin-33, to elucidate the phenotype of recipient cells in the context of uptake-independent ligand-receptor interactions.

DECLARATIONS

Authors' contributions

The author contributed solely to the article

Availability of data and materials

Not applicable.

Financial support and sponsorship

None.

Conflicts of interest

Pfaffl MW is an Editorial Board member of the journal *Extracellular Vesicles and Circulating Nucleic Acids*.

Ethical approval and consent to participate

Not applicable.

Consent for publication

Not applicable.

Copyright

© The Author(s) 2024.

REFERENCES

1. Welsh JA, Goberdhan DCI, O'Driscoll L, et al; MISEV Consortium. Minimal information for studies of extracellular vesicles (MISEV2023): from basic to advanced approaches. *J Extracell Vesicles* 2024;13:e12404. DOI PubMed PMC
2. Corrigan L, Redhai S, Leiblich A, et al. BMP-regulated exosomes from *Drosophila* male reproductive glands reprogram female behavior. *J Cell Biol* 2014;206:671-88. DOI PubMed PMC

3. Fan SJ, Kroeger B, Marie PP, et al. Glutamine deprivation alters the origin and function of cancer cell exosomes. *EMBO J* 2020;39:e103009. [DOI](#) [PubMed](#) [PMC](#)
4. Verweij FJ, Balaj L, Boulanger CM, et al. The power of imaging to understand extracellular vesicle biology in vivo. *Nat Methods* 2021;18:1013-26. [DOI](#) [PubMed](#) [PMC](#)
5. Kur IM, Prouvot PH, Fu T, et al. Neuronal activity triggers uptake of hematopoietic extracellular vesicles in vivo. *PLoS Biol* 2020;18:e3000643. [DOI](#) [PubMed](#) [PMC](#)
6. Zomer A, Maynard C, Verweij FJ, et al. In Vivo imaging reveals extracellular vesicle-mediated phenocopying of metastatic behavior. *Cell* 2015;161:1046-57. [DOI](#) [PubMed](#) [PMC](#)
7. Zomer A, Steenbeek SC, Maynard C, van Rheeën J. Studying extracellular vesicle transfer by a Cre-loxP method. *Nat Protoc* 2016;11:87-101. [DOI](#) [PubMed](#)
8. Fordjour FK, Abuelreich S, Hong X, et al. Exomap1 mouse: a transgenic model for *in vivo* studies of exosome biology. *bioRxiv* 2023:Online ahead of print. [DOI](#) [PubMed](#) [PMC](#)
9. Alashkar Alhamwe B, Ponath V, Alhamdan F, et al. BAG6 restricts pancreatic cancer progression by suppressing the release of IL33-presenting extracellular vesicles and the activation of mast cells. *Cell Mol Immunol* 2024;21:918-31. [DOI](#) [PubMed](#) [PMC](#)
10. Pham TND, Shields MA, Spaulding C, et al. Preclinical models of pancreatic ductal adenocarcinoma and their utility in immunotherapy studies. *Cancers* 2021;13:440. [DOI](#) [PubMed](#) [PMC](#)
11. Wang Y, Zhang Y, Yang J, et al. Genomic sequencing of key genes in mouse pancreatic cancer cells. *Curr Mol Med* 2012;12:331-41. [DOI](#) [PubMed](#) [PMC](#)
12. Schuldner M, Dörsam B, Shatnyeva O, et al. Exosome-dependent immune surveillance at the metastatic niche requires BAG6 and CBP/p300-dependent acetylation of p53. *Theranostics* 2019;9:6047-62. [DOI](#) [PubMed](#) [PMC](#)